

NLRP7 Antibody, HRP conjugated



PACO57353

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Inhibits CASP1/caspase-1-dependent IL1B secretion.

Gene ID:

NLRP7

Uniprot

Q8WX94

Synonyms:

NACHT, LRR and PYD domains-containing protein 7 (Nucleotide-binding oligomerization domain protein 12) (PYRIN-containing APAF1-like protein 3), NLRP7, NALP7 NOD12 PYPAF3

Immunogen:

Recombinant Human NACHT, LRR and PYD domains-containing protein 7 protein (89-187AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A