

# NLGN4X Antibody



PACO57364

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IF:1:200-1:500

**Protein Background:**

Putative neuronal cell surface protein involved in cell-cell-interactions.

**Gene ID:**

NLGN4X

**Uniprot**

Q8N0W4

**Synonyms:**

Neurologin-4, X-linked (Neurologin X) (HNLX), NLGN4X, KIAA1260 NLGN4

**Immunogen:**


Recombinant Human Neurologin-4, X-linked protein (383-494AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Immunofluorescence staining of HeLa cells with PACO57364 at 1:200, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).