

# NLGN4X Antibody, HRP conjugated



PACO57365

---

## Product Information

<b>Size:</b>	<b>Protein Background:</b>
50ug	Putative neuronal cell surface protein involved in cell-cell-interactions.
<b>Reactivity:</b>	<b>Gene ID:</b>
Human	NLGN4X
<b>Source:</b>	<b>Uniprot</b>
Rabbit	Q8N0W4
<b>Isotype:</b>	<b>Synonyms:</b>
IgG	Neurologin-4, X-linked (Neurologin X) (HNLX), NLGN4X, KIAA1260 NLGN4
<b>Applications:</b>	<b>Immunogen:</b>
ELISA	Recombinant Human Neurologin-4, X-linked protein (383-494AA).
<b>Recommended dilutions:</b>	<b>Storage:</b>
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A