

# DDX17 Antibody, HRP conjugated



PACO57369

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

As an RNA helicase, unwinds RNA and alters RNA structures through ATP binding and hydrolysis. Involved in multiple cellular processes, including pre-mRNA splicing, alternative splicing, ribosomal RNA processing and miRNA processing, as well as transcription regulation. Regulates the alternative splicing of exons exhibiting specific features. For instance, promotes the inclusion of AC-rich alternative exons in CD44 transcripts. This function requires the RNA helicase activity.

**Gene ID:**

DDX17

**Uniprot**

Q92841

**Synonyms:**

Probable ATP-dependent RNA helicase DDX17 (EC 3.6.4.13) (DEAD box protein 17) (DEAD box protein p72) (DEAD box protein p82) (RNA-dependent helicase p72), DDX17

**Immunogen:**

Recombinant Human Probable ATP-dependent RNA helicase DDX17 protein (591-692AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A