

THEMIS Antibody, FITC conjugated



PACO57386

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a central role in late thymocyte development by controlling both positive and negative T-cell selection. Required to sustain and/or integrate signals required for proper lineage commitment and maturation of T-cells. Regulates T-cell development through T-cell antigen receptor (TCR) signaling and in particular through the regulation of calcium influx and phosphorylation of Erk.

Gene ID:

THEMIS

Uniprot

Q8N1K5

Synonyms:

Protein THEMIS (Thymocyte-expressed molecule involved in selection), THEMIS, C6orf190 C6orf207

Immunogen:

Recombinant Human Protein THEMIS protein (547-641AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A