

AQP10 Antibody, Biotin conjugated



PACO57403

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Water channel required to promote glycerol permeability and water transport across cell membranes. May contribute to water transport in the upper portion of small intestine. Isoform 2 is not permeable to urea and glycerol.

Gene ID:

AQP10

Uniprot

Q96PS8

Synonyms:

Aquaporin-10 (AQP-10) (Aquaglyceroporin-10) (Small intestine aquaporin), AQP10

Immunogen:

Recombinant Human Aquaporin-10 protein (209-301AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A