

SLC26A6 Antibody, HRP conjugated



PACO57433

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

SLC26A6

Uniprot

Q9BXS9

Synonyms:

Solute carrier family 26 member 6 (Anion exchange transporter) (Pendrin-like protein 1) (Pendrin-L1), SLC26A6

Immunogen:

Recombinant Human Solute carrier family 26 member 6 protein (554-668AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A