

SCUBE2 Antibody, FITC conjugated



PACO57450

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cell surface, extracellular region, extracellular space, cell differentiation, signal transduction.

Gene ID:

SCUBE2

Uniprot

Q9NQ36

Synonyms:

Signal peptide, CUB and EGF-like domain-containing protein 2 (Protein CEGP1) (Scube/You), SCUBE2, CEGP1

Immunogen:

Recombinant Human Signal peptide, CUB and EGF-like domain-containing protein 2 protein (534-703AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A