

SEL1L Antibody, Biotin conjugated



PACO57455

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in the endoplasmic reticulum quality control (ERQC) system also called ER-associated degradation (ERAD) involved in ubiquitin-dependent degradation of misfolded endoplasmic reticulum proteins. Enhances SYVN1 stability. Plays a role in LPL maturation and secretion. Required for normal differentiation of the pancreas epithelium, and for normal exocrine function and survival of pancreatic cells. May play a role in Notch signaling.

Gene ID:

SEL1L

Uniprot

Q9UBV2

Synonyms:

Protein sel-1 homolog 1 (Suppressor of lin-12-like protein 1) (Sel-1L), SEL1L, TSA305

Immunogen:

Recombinant Human Protein sel-1 homolog 1 protein (22-181AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A