

DDX28 Antibody, HRP conjugated



PACO57501

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays an essential role in facilitating the proper assembly of the mitochondrial large ribosomal subunit and its helicase activity is essential for this function. May be involved in RNA processing or transport. Has RNA and Mg(2+)-dependent ATPase activity.

Gene ID:

DDX28

Uniprot

Q9NUL7

Synonyms:

Probable ATP-dependent RNA helicase DDX28 (EC 3.6.4.13) (Mitochondrial DEAD box protein 28), DDX28, MDDX28

Immunogen:

Recombinant Human Probable ATP-dependent RNA helicase DDX28 protein (71-206AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A