

# ONECUT1 Antibody, HRP conjugated



PACO57509

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transcriptional activator. Binds the consensus sequence 5'-DHWATTGAYTWWD-3' on a variety of gene promoters such as those of HNF3B and TTR. Important for liver genes transcription.

**Gene ID:**

ONECUT1

**Uniprot**

Q9UBC0

**Synonyms:**

Hepatocyte nuclear factor 6 (HNF-6) (One cut domain family member 1) (One cut homeobox 1), ONECUT1, HNF6 HNF6A

**Immunogen:**

Recombinant Human Hepatocyte nuclear factor 6 protein (192-292AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A