

FANCL Antibody, FITC conjugated



PACO57530

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Ubiquitin ligase protein that mediates monoubiquitination of FANCD2, a key step in the DNA damage pathway. Also mediates monoubiquitination of FANCI. May stimulate the ubiquitin release from UBE2W. May be required for proper primordial germ cell proliferation in the embryonic stage, whereas it is probably not needed for spermatogonial proliferation after birth.

Gene ID:

FANCL

Uniprot

Q9NW38

Synonyms:

E3 ubiquitin-protein ligase FANCL (EC 2.3.2.27) (Fanconi anemia group L protein) (Fanconi anemia-associated polypeptide of 43 kDa) (FAAP43) (RING-type E3 ubiquitin transferase FANCL), FANCL, PHF9

Immunogen:

Recombinant Human E3 ubiquitin-protein ligase FANCL protein (138-248AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A