

DMGDH Antibody, HRP conjugated



PACO57533

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the demethylation of N, N-dimethylglycine to sarcosine. Also has activity with sarcosine in vitro.

Gene ID:

DMGDH

Uniprot

Q9UI17

Synonyms:

Dimethylglycine dehydrogenase, mitochondrial (EC 1.5.8.4) (ME2GLYDH), DMGDH

Immunogen:

Recombinant Human Dimethylglycine dehydrogenase, mitochondrial protein (429-524AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A