

# DMGDH Antibody, Biotin conjugated



PACO57535

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Catalyzes the demethylation of N, N-dimethylglycine to sarcosine. Also has activity with sarcosine in vitro.

**Gene ID:**

DMGDH

**Uniprot**

Q9UI17

**Synonyms:**

Dimethylglycine dehydrogenase, mitochondrial (EC 1.5.8.4) (ME2GLYDH), DMGDH

**Immunogen:**

Recombinant Human Dimethylglycine dehydrogenase, mitochondrial protein (429-524AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A