

USP24 Antibody, HRP conjugated



PACO57569

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Protease that can remove conjugated ubiquitin from target proteins and polyubiquitin chains. Deubiquitinates DDB2, preventing its proteasomal degradation.

Gene ID:

USP24

Uniprot

Q9UPU5

Synonyms:

Ubiquitin carboxyl-terminal hydrolase 24 (EC 3.4.19.12) (Deubiquitinating enzyme 24) (Ubiquitin thioesterase 24) (Ubiquitin-specific-processing protease 24), USP24, KIAA1057

Immunogen:

Recombinant Human Ubiquitin carboxyl-terminal hydrolase 24 protein (2509-2620AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A