

HAO1 Antibody, FITC conjugated



PACO57578

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Has 2-hydroxyacid, oxidase activity. Most active on the 2-carbon substrate glycolate, but is also active on 2-hydroxy fatty acid, , with high activity towards 2-hydroxy palmitate and 2-hydroxy octanoate.

Gene ID:

HAO1

Uniprot

Q9UJM8

Synonyms:

Hydroxyacid, oxidase 1 (HAOX1) (EC 1.1.3.15) (Glycolate oxidase) (GOX), HAO1, GOX1
HAOX1

Immunogen:

Recombinant Human Hydroxyacid, oxidase 1 protein (161-273AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A