

CDH9 Antibody, FITC conjugated



PACO57586

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cadherins are calcium-dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

Gene ID:

CDH9

Uniprot

Q9ULB4

Synonyms:

Cadherin-9, CDH9

Immunogen:

Recombinant Human Cadherin-9 protein (131-314AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A