

SOX21 Antibody, FITC conjugated



PACO57606

Product Information

Size:	Protein Background:
50ug	May play a role as an activator of transcription of OPRM1.
Reactivity:	Gene ID:
Human	SOX21
Source:	Uniprot
Rabbit	Q9Y651
Isotype:	Synonyms:
IgG	Transcription factor SOX-21 (SOX-A), SOX21, SOX25 SOXA
Applications:	Immunogen:
ELISA	Recombinant Human Transcription factor SOX-21 protein (25-126AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A