

# RAD54L2 Antibody, HRP conjugated



PACO57609

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

DNA helicase that modulates androgen receptor (AR)-dependent transactivation in a promoter-dependent manner. Not able to remodel mononucleosomes in vitro.

**Gene ID:**

RAD54L2

**Uniprot**

Q9Y4B4

**Synonyms:**

Helicase ARIP4 (EC 3.6.4.12) (Androgen receptor-interacting protein 4) (RAD54-like protein 2), RAD54L2, ARIP4 KIAA0809

**Immunogen:**

Recombinant Human Helicase ARIP4 protein (1122-1296AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A