

RAD54L2 Antibody, Biotin conjugated



PACO57611

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

DNA helicase that modulates androgen receptor (AR)-dependent transactivation in a promoter-dependent manner. Not able to remodel mononucleosomes in vitro.

Gene ID:

RAD54L2

Uniprot

Q9Y4B4

Synonyms:

Helicase ARIP4 (EC 3.6.4.12) (Androgen receptor-interacting protein 4) (RAD54-like protein 2), RAD54L2, ARIP4 KIAA0809

Immunogen:

Recombinant Human Helicase ARIP4 protein (1122-1296AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A