

APOBEC1 Antibody, HRP conjugated



PACO57661

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalytic component of the apolipoprotein B mRNA editing enzyme complex which is responsible for the posttranscriptional editing of a CAA codon for Gln to a UAA codon for stop in the APOB mRNA. Also involved in CGA (Arg) to UGA (Stop) editing in the NF1 mRNA. May also play a role in the epigenetic regulation of gene expression by participating in DNA demethylation.

Gene ID:

APOBEC1

Uniprot

P41238

Synonyms:

C->U-editing enzyme APOBEC-1 (EC 3.5.4.36) (Apolipoprotein B mRNA-editing enzyme 1) (HEPR) (mRNA(cytosine(6666)) deaminase 1), APOBEC1

Immunogen:

Recombinant Human C->U-editing enzyme APOBEC-1 protein (1-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A