

BTD Antibody, FITC conjugated



PACO57678

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalytic release of biotin from biocytin, the product of biotin-dependent carboxylases degradation.

Gene ID:

BTD

Uniprot

P43251

Synonyms:

Biotinidase (Biotinase) (EC 3.5.1.12), BTD

Immunogen:

Recombinant Human Biotinidase protein (322-397AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A