

# CRLF1 Antibody, FITC conjugated



PACO57722

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Cytokine receptor subunit, possibly playing a regulatory role in the immune system and during fetal development. May be involved in nervous system development.

**Gene ID:**

CRLF1

**Uniprot**

O75462

**Synonyms:**

Cytokine receptor-like factor 1 (Cytokine-like factor 1) (CLF-1) (ZcytoR5), CRLF1

**Immunogen:**

Recombinant Human Cytokine receptor-like factor 1 protein (318-422AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A