

DNPEP Antibody, HRP conjugated



PACO57729

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Aminopeptidase with specificity towards an acid, c amino acid, at the N-terminus. Likely to play an important role in intracellular protein and peptide metabolism.

Gene ID:

DNPEP

Uniprot

Q9ULA0

Synonyms:

Aspartyl aminopeptidase (EC 3.4.11.21), DNPEP, ASPEP DAP

Immunogen:

Recombinant Human Aspartyl aminopeptidase protein (279-386AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A