

# DNPEP Antibody, Biotin conjugated



PACO57731

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Aminopeptidase with specificity towards an acid, c amino acid, at the N-terminus. Likely to play an important role in intracellular protein and peptide metabolism.

**Gene ID:**

DNPEP

**Uniprot**

Q9ULA0

**Synonyms:**

Aspartyl aminopeptidase (EC 3.4.11.21), DNPEP, ASPEP DAP

**Immunogen:**

Recombinant Human Aspartyl aminopeptidase protein (279-386AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A