

DPPA5 Antibody, FITC conjugated



PACO57734

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in the maintenance of embryonic stem (ES) cell pluripotency. Dispensable for self-renewal of pluripotent ES cells and establishment of germ cells. Associates with specific target mRNAs.

Gene ID:

DPPA5

Uniprot

A6NC42

Synonyms:

Developmental pluripotency-associated 5 protein (hDPPA5) (Embryonal stem cell-specific gene 1 protein) (ESG-1), DPPA5, ESG1

Immunogen:

Recombinant Human Developmental pluripotency-associated 5 protein (1-116AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A