

# FOXL2 Antibody, Biotin conjugated



PACO57771

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transcriptional regulator. Critical factor essential for ovary differentiation and maintenance, and repression of the genetic program for somatic testis determination. Prevents trans-differentiation of ovary to testis through transcriptional repression of the Sertoli cell-promoting gene SOX9. Has apoptotic activity in ovarian cells. Suppresses ESR1-mediated transcription of PTGS2/COX2 stimulated by tamoxifen. Is a regulator of CYP19 expression. Participates in SMAD3-dependent transcription of FST via the intronic SMAD-binding element. Is a transcriptional repressor of STAR. Activates SIRT1 transcription under cellular stress conditions. Activates transcription of OSR2.

**Gene ID:**

FOXL2

**Uniprot**

P58012

**Synonyms:**

Forkhead box protein L2, FOXL2

**Immunogen:**

Recombinant Human Forkhead box protein L2 protein (116-216AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A