

GBX2 Antibody, HRP conjugated



PACO57773

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May act as a transcription factor for cell pluripotency and differentiation in the embryo.

Gene ID:

GBX2

Uniprot

P52951

Synonyms:

Homeobox protein GBX-2 (Gastrulation and brain-specific homeobox protein 2), GBX2

Immunogen:

Recombinant Human Homeobox protein GBX-2 protein (65-248AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A