

GTPBP1 Antibody, FITC conjugated



PACO57798

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Promotes degradation of target mRNA species. Plays a role in the regulation of circadian mRNA stability. Binds GTP and has GTPase activity.

Gene ID:

GTPBP1

Uniprot

O00178

Synonyms:

GTP-binding protein 1 (G-protein 1) (GP-1) (GP1), GTPBP1

Immunogen:

Recombinant Human GTP-binding protein 1 protein (562-669AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A