

# INPP5B Antibody, HRP conjugated



PACO57821

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Hydrolyzes phosphatidylinositol 4,5-bisphosphate (PtdIns(4,5)P<sub>2</sub>) and the signaling molecule phosphatidylinositol 1,4,5-trisphosphate (PtdIns(1,4,5)P<sub>3</sub>), and thereby modulates cellular signaling events.

**Gene ID:**

INPP5B

**Uniprot**

P32019

**Synonyms:**

Type II inositol 1,4,5-trisphosphate 5-phosphatase (EC 3.1.3.36) (75 kDa inositol polyphosphate-5-phosphatase) (Phosphoinositide 5-phosphatase) (5PTase), INPP5B, OCRL2

**Immunogen:**

Recombinant Human Type II inositol 1,4,5-trisphosphate 5-phosphatase protein (78-178AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A