

KRT1 Antibody, Biotin conjugated



PACO57835

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May regulate the activity of kinases such as PKC and SRC via binding to integrin beta-1 (ITB1) and the receptor of activated protein C kinase 1 (RACK1). In complex with C1QBP is a high affinity receptor for kininogen-1/HMWK.

Gene ID:

KRT1

Uniprot

P04264

Synonyms:

Keratin, type II cytoskeletal 1 (67 kDa cyokeratin) (Cytokeratin-1) (CK-1) (Hair alpha protein) (Keratin-1) (K1) (Type-II keratin Kb1), KRT1, KRTA

Immunogen:

Recombinant Human Keratin, type II cytoskeletal 1 protein (267-438AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A