

# LIFR Antibody, Biotin conjugated



PACO57843

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Signal-transducing molecule. May have a common pathway with IL6ST. The soluble form inhibits the biological activity of LIF by blocking its binding to receptors on target cells.

**Gene ID:**

LIFR

**Uniprot**

P42702

**Synonyms:**

Leukemia inhibitory factor receptor (LIF receptor) (LIF-R) (CD antigen CD118), LIFR

**Immunogen:**

Recombinant Human Leukemia inhibitory factor receptor protein (915-1086AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A