

# LSAMP Antibody, FITC conjugated



PACO57850

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Mediates selective neuronal growth and axon targeting. Contributes to the guidance of developing axons and remodeling of mature circuits in the limbic system. Essential for normal growth of the hippocampal mossy fiber projection.

**Gene ID:**

LSAMP

**Uniprot**

Q13449

**Synonyms:**

Limbic system-associated membrane protein (LSAMP) (IGLON family member 3), LSAMP, IGLON3 LAMP

**Immunogen:**

Recombinant Human Limbic system-associated membrane protein (168-287AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A