

LSAMP Antibody, Biotin conjugated



PACO57851

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Mediates selective neuronal growth and axon targeting. Contributes to the guidance of developing axons and remodeling of mature circuits in the limbic system. Essential for normal growth of the hippocampal mossy fiber projection.

Gene ID:

LSAMP

Uniprot

Q13449

Synonyms:

Limbic system-associated membrane protein (LSAMP) (IGLON family member 3), LSAMP, IGLON3 LAMP

Immunogen:

Recombinant Human Limbic system-associated membrane protein (168-287AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A