

MC3R Antibody, FITC conjugated



PACO57854

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Receptor for MSH (alpha, beta and gamma) and ACTH. This receptor is mediated by G proteins which activate adenylate cyclase. Required for expression of anticipatory patterns of activity and wakefulness during periods of limited nutrient availability and for the normal regulation of circadian clock activity in the brain.

Gene ID:

MC3R

Uniprot

P41968

Synonyms:

Melanocortin receptor 3 (MC3-R), MC3R

Immunogen:

Recombinant Human Melanocortin receptor 3 protein (38-74AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A