

# MVD Antibody, HRP conjugated



PACO57869

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Performs the first committed step in the biosynthesis of isoprenes.

**Gene ID:**

MVD

**Uniprot**

P53602

**Synonyms:**

Diphosphomevalonate decarboxylase (EC 4.1.1.33) (Mevalonate (diphospho)decarboxylase) (MDDase) (Mevalonate pyrophosphate decarboxylase), MVD, MPD

**Immunogen:**

Recombinant Human Diphosphomevalonate decarboxylase protein (69-221AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A