

MYH11 Antibody, Biotin conjugated



PACO57879

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Muscle contraction.

Gene ID:

MYH11

Uniprot

P35749

Synonyms:

Myosin-11 (Myosin heavy chain 11) (Myosin heavy chain, smooth muscle isoform) (SMMHC), MYH11, KIAA0866

Immunogen:

Recombinant Human Myosin-11 protein (1655-1797AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A