

POLD2 Antibody, FITC conjugated



PACO57926

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	POLD2
Source:	Uniprot
Rabbit	P49005
Isotype:	Synonyms:
IgG	DNA polymerase delta subunit 2 (DNA polymerase delta subunit p50), POLD2
Applications:	Immunogen:
ELISA	Recombinant Human DNA polymerase delta subunit 2 protein (245-430AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A