

# PTPN3 Antibody, FITC conjugated



PACO57934

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May act at junctions between the membrane and the cytoskeleton. Possesses tyrosine phosphatase activity.

**Gene ID:**

PTPN3

**Uniprot**

P26045

**Synonyms:**

Tyrosine-protein phosphatase non-receptor type 3 (EC 3.1.3.48) (Protein-tyrosine phosphatase H1) (PTP-H1), PTPN3, PTPH1

**Immunogen:**

Recombinant Human Tyrosine-protein phosphatase non-receptor type 3 protein (365-505AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A