

# RFXANK Antibody, FITC conjugated



PACO57946

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Activates transcription from class II MHC promoters. Activation requires the activity of the MHC class II transactivator/CIITA. May regulate other genes in the cell. RFX binds the X1 box of MHC-II promoters. May also potentiate the activation of RAF1.

**Gene ID:**

RFXANK

**Uniprot**

O14593

**Synonyms:**

DNA-binding protein RFXANK (Ankyrin repeat family A protein 1) (Regulatory factor X subunit B) (RFX-B) (Regulatory factor X-associated ankyrin-containing protein), RFXANK, ANKRA1 RFXB

**Immunogen:**

Recombinant Human DNA-binding protein RFXANK protein (3-84AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A