

RFXANK Antibody, Biotin conjugated



PACO57947

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Activates transcription from class II MHC promoters. Activation requires the activity of the MHC class II transactivator/CIITA. May regulate other genes in the cell. RFX binds the X1 box of MHC-II promoters. May also potentiate the activation of RAF1.

Gene ID:

RFXANK

Uniprot

O14593

Synonyms:

DNA-binding protein RFXANK (Ankyrin repeat family A protein 1) (Regulatory factor X subunit B) (RFX-B) (Regulatory factor X-associated ankyrin-containing protein), RFXANK, ANKRA1 RFXB

Immunogen:

Recombinant Human DNA-binding protein RFXANK protein (3-84AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A