

# SOX8 Antibody, HRP conjugated



PACO57985

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role in central nervous system, limb and facial development. May be involved in male sex determination. Binds the consensus motif 5'-[AT][AT]CAA[AT]G-3'.

**Gene ID:**

SOX8

**Uniprot**

P57073

**Synonyms:**

Transcription factor SOX-8, SOX8

**Immunogen:**

Recombinant Human Transcription factor SOX-8 protein (292-386AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A