

# SRGAP3 Antibody, FITC conjugated



PACO57994

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

GTPase-activating protein for RAC1 and perhaps Cdc42, but not for RhoA small GTPase. May attenuate RAC1 signaling in neurons.

**Gene ID:**

SRGAP3

**Uniprot**

O43295

**Synonyms:**

SLIT-ROBO Rho GTPase-activating protein 3 (srGAP3) (Mental disorder-associated GAP) (Rho GTPase-activating protein 14) (WAVE-associated Rac GTPase-activating protein) (WRP), SRGAP3, ARHGAP14 KIAA0411 KIAA1156 MEGAP SRGAP2

**Immunogen:**

Recombinant Human SLIT-ROBO Rho GTPase-activating protein 3 protein (709-955AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A