

# TRAPPC6A Antibody, FITC conjugated



PACO58010

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play a role in vesicular transport during the biogenesis of melanosomes.

**Gene ID:**

TRAPPC6A

**Uniprot**

O75865

**Synonyms:**

Trafficking protein particle complex subunit 6A (TRAPP complex subunit 6A), TRAPPC6A

**Immunogen:**

Recombinant Human Trafficking protein particle complex subunit 6A protein (1-173AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A