

UBE2G1 Antibody, HRP conjugated



PACO58017

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-48'-, as well as 'Lys-63'-linked polyubiquitination. May be involved in degradation of muscle-specific proteins. Mediates polyubiquitination of CYP3A4.

Gene ID:

UBE2G1

Uniprot

P62253

Synonyms:

Ubiquitin-conjugating enzyme E2 G1 (EC 2.3.2.23) (E2 ubiquitin-conjugating enzyme G1) (E217K) (UBC7) (Ubiquitin carrier protein G1) (Ubiquitin-protein ligase G1) [Cleaved into: Ubiquitin-conjugating enzyme E2 G1, N-terminally processed], UBE2G1, UBE2G

Immunogen:

Recombinant Human Ubiquitin-conjugating enzyme E2 G1 protein (2-170AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A