

ZFX Antibody, Biotin conjugated



PACO58051

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable transcriptional activator.

Gene ID:

ZFX

Uniprot

P17010

Synonyms:

Zinc finger X-chromosomal protein, ZFX

Immunogen:

Recombinant Human Zinc finger X-chromosomal protein (711-805AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A