

ADAM8 Antibody, FITC conjugated



PACO58058

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Possible involvement in extravasation of leukocytes.

Gene ID:

ADAM8

Uniprot

P78325

Synonyms:

Disintegrin and metalloproteinase domain-containing protein 8 (ADAM 8) (EC 3.4.24)
(Cell surface antigen MS2) (CD antigen CD156a), ADAM8, MS2

Immunogen:

Recombinant Human Disintegrin and metalloproteinase domain-containing protein 8 protein (677-824AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A