

DPY19L1 Antibody



PACO58076

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Probable C-mannosyltransferase that mediates C-mannosylation of tryptophan residues on target proteins.

Gene ID:

DPY19L1

Uniprot

Q2PZI1

Synonyms:

Probable C-mannosyltransferase DPY19L1 (EC 2.4.1) (Dpy-19-like protein 1) (Protein dpy-19 homolog 1), DPY19L1, GA0500 KIAA0877


Immunogen:

Recombinant Human Probable C-mannosyltransferase DPY19L1 protein (513-665AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescence staining of A549 cells with PACO58076 at 1:166, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).