

# CUX1 Antibody, HRP conjugated



PACO58121

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May be involved in intra-Golgi retrograde transport.

**Gene ID:**

CUX1

**Uniprot**

Q13948

**Synonyms:**

Protein CASP, CUX1, CUTL1

**Immunogen:**

Recombinant Human Protein CASP protein (50-210AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A