

FASTKD5 Antibody, Biotin conjugated



PACO58219

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays an important role in the processing of non-canonical mitochondrial mRNA precursors.

Gene ID:

FASTKD5

Uniprot

Q7L8L6

Synonyms:

FAST kinase domain-containing protein 5, mitochondrial, FASTKD5, KIAA1792

Immunogen:

Recombinant Human FAST kinase domain-containing protein 5, mitochondrial protein (516-664AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A